

ADHESIVE TAPES

A full range


Amcor Flexibles SPS offers a full range of neutral and indicator adhesive tapes, adapted to steam sterilisation and ethylene oxide gas processes.

The type of adhesive used has been designed specially to prevent any migration during the sterilisation process. Amcor Flexibles SPS adhesive tapes do not contain any lead and pose no risk of toxicity for the packaged devices or the user personnel. The tapes contain no latex in their composition.

The tapes come in several widths and are available with two types of adhesiveness, depending on the kind of material on which they are to be applied:

- White for use on "Crepe" material
- Blue for use on "Non-woven" material

All our tapes are available either in neutral form or with flow indicators for sterilisation.



FOR CREPE PAPER (beige)

Neutral version:

Steam autoclave sterilisation
Width: 18mm / 24mm / 48mm

Steam indicator version:

Tell-tale indicator
Width: 18mm / 24mm

FOR NON-WOVEN MATERIAL (blue)

For maximum adhesion, specially adapted to "non-woven" material.

Neutral version:

Width: 24mm

Steam indicator version:

Tell-tale indicator
Width: 18mm / 24mm

COMPLIANCE

EN ISO 11 140-1 (Indicator)

ETHYLENE OXIDE INDICATOR TAPE (beige)

Adhesive mass adapted to the EO gas sterilisation process
Width: 18mm