



handicare

Making everyday life easier

# MiniLift 125 / 160 / 200

A totally unique and very comfortable sit-to-stand experience. Together with the MiniLift, the ThoraxSling with seat support provides a very gentle and pleasant sit-to-stand experience.

The movement of the lift arm is designed so that the pulling motion ends comfortably in the lower back, which prevents the vest from riding up under the arms.

Available in lifting weights of 125kg / 160kg / 200kg and with manual or electrical based - width adjustments (models with max lifting weight of 160kg only).



## ThoraxSling with seat support

- Specifically adapted sit-to-stand vest that is used in combination with the MiniLift
- Easy for the caregiver to apply
- Very comfortable and gives the user an added sense of security
- Facilitates toileting
- Does not ride up under the arms
- Alternative lift loops enable an upright posture, regardless the user's height



Thorax Sling with seat support

|                         |          |          |           |            |
|-------------------------|----------|----------|-----------|------------|
| <i>Code</i>             | 45600004 | 45600006 | 45600007  | 45600008   |
| <i>Size</i>             | S        | M        | L         | XL         |
| <i>Rec. user weight</i> | 30-60 kg | 50-80 kg | 70-110 kg | 100-150 kg |
| <i>Lift vest width</i>  | 94 cm    | 104 cm   | 124 cm    | 146 cm     |
| <i>SWL</i>              | 300kg    | 300kg    | 300kg     | 300kg      |

For all order enquiries please email [admin@medix21.co.nz](mailto:admin@medix21.co.nz) or Tel 04 233 1627

Exclusive to:  
Medix 21  
29 Jarden Mile  
Ngauranga 6035  
Wellington

