

**HARTMANN**



# Continence Management Catalogue

2023

**MoliCare®**





# MoliCare®

## Protection. Every day.

HARTMANN has been providing healthcare products for over 200 years and incontinence products for over 50 years. We take a holistic approach to healthcare, placing emphasis on the tangible benefits for patients and those who care for them, aiming to enhance their well-being.

## The MoliCare Difference

MoliCare is dedicated to meeting the most stringent standards of quality, dependability, and effectiveness by utilising our advanced technology to develop innovative solutions for managing incontinence. The aim is to address the needs of individuals living with incontinence, enhance the daily duties of skilled caregivers, and satisfy the expectations of healthcare professionals.

Incontinence can disrupt the skin's natural barrier and acid protection mantle, leaving it more susceptible to irritation and damage from urine and faeces.<sup>1</sup> In order to prevent Incontinence Associated Dermatitis (IAD), which can cause skin inflammation, blisters, itchiness or skin erosion<sup>2-4</sup>, it is important to take measures to protect the skin from the harmful effects of incontinence.

MoliCare provides absorbent products that are specifically designed to keep the skin dry and maintain its acid protection layer. This can help alleviate the negative consequences of stress factors associated with incontinence and prevent the development of IAD. This is in line with expert recommendations, which emphasises the importance of preventing IAD as a means of improving quality of life for individuals living with incontinence.<sup>5</sup>



pH of Skin on Incontinence Briefs<sup>6</sup>



### MoliCare Absorbent Products:


- ✓ Significantly lower skin pH values compared to standard briefs, when skin is exposed to alkaline solution.<sup>6,7</sup>
- ✓ Lower the prevalence of erythema<sup>6,7</sup> which is skin redness, an early sign of skin irritation.

Graph data from Bliss Study Absorbent Briefs Containing Curly Fiber Lower (Acidify) Skin pH Reducing Risk for Incontinence Associated Skin Damage in Older Nursing Home Residents.<sup>6</sup> ANOVA F=8.83 for thigh and F=10.12 for forearm pH differences, p < .001 for both; in post-hoc comparisons, skin pH values after wet brief using MoliCare are lower than all other conditions (p=.001 for all).<sup>6</sup>




### ✓ MoliCare controls odour and promotes a healthy skin environment.

- Buffers urine to a skin-friendly pH
- Maintains a skin-healthy pH balance
- Preserves the acid protection mantle
- Facilitates rapid fluid distribution within the brief
- Keeps the skin feeling dry



View MoliCare's optimal skin compatibility video to understand how it works.

Scan the code or visit:  
[lead.me/ph-balance](https://lead.me/ph-balance)



1. Shigeta, Yoshie, et al. "Exploring the relationship between skin property and absorbent pad environment." Journal of Clinical Nursing 18.11 (2009): 1607-1616. 2. Gray, Mikel. "Optimal management of incontinence-associated dermatitis in the elderly." American journal of clinical dermatology 11.3 (2010): 201-210. 3. Gray, Mikel, et al. "Incontinence-associated dermatitis: a consensus." Journal of Wound Ostomy & Continence Nursing 34.1 (2007): 45-54. 4. Junkin, Joan, and Joan Lerner Selekof. "Beyond 'diaper rash': Incontinence-associated dermatitis Does it have you seeing RED?" Nursing2020 38.11 (2008): 56hn1-56hn10. Warner RR, Stone KJ, Boissy YL. Hydration disrupts human stratum corneum ultrastructure. J Invest Dermatol. 2003;120(2):275-84. 5. Gray, Mikel, et al. "Incontinence-associated dermatitis: a comprehensive review and update." Journal of Wound Ostomy & Continence Nursing 39.1 (2012): 61-74. 6. Bliss et al., An economic evaluation of four skin damage prevention regimens in nursing home residents with incontinence; economics of skin damage prevention; J Wound Ostomy Continence Nurse. 2007. 7. Data on file. The protective effect of the new diaper design on the epidermal barrier. 2022.






# Health Economics of MoliCare

**Incidence**  
 In an acute care study of 5,432 patients<sup>1</sup>:  
 • 45.7% of patients had a prevalence of IAD  
 • 73% of IAD patients developed IAD following admission  
 IAD accompanied by unrelieved pressure results in a 5-fold increase in pressure injury risk<sup>2</sup>

**Cost - IAD<sup>1</sup>**  
 The average cost to treat one mild to moderate IAD episode in one resident is \$220.  
**IAD occurrence can be up to 52%, resulting in costs of \$3,294 per month.**  
 (This example assumes 20% in a 100 bed long term care facility with 75% incontinence at April 2018 rate<sup>1</sup>)




1. Gray,M., Giuliano, K.K., JWOCN, Incontinence Associated Dermatitis, Characteristics and Relationship to Pressure Injury A Multisite Epidemiologic Analysis, January/February 2018. 2. Melter, Pressure Ulcer Assessment and Treatment; Wild Iris Medical Education 2014; 17. 3. Bliss et al., An economic evaluation of four skin damage prevention regimens in nursing home residents with incontinence; economics of skin damage prevention; J Wound Ostomy Continence Nurse. 2007;34(2): 143-152.

# The MoliCare Skintegrity System

MoliCare offers a comprehensive approach to incontinence management by combining MoliCare absorbent products with MoliCare Skin care products. This combination, known as the MoliCare Skintegrity System, helps to maintain a healthy skin pH balance and help prevent the occurrence of Incontinence Associated Dermatitis (IAD).

MoliCare absorbent products work in conjunction with MoliCare Skin care products to maintain the acid protection mantle and provide a skin-neutral pH balance.<sup>1</sup>  
 Our goal is to provide innovative solutions that improve the management of incontinence, allowing individuals to live with greater comfort and dignity.



1. Except MoliCare® Premium Men Pad 3 Drops, MoliCare® Skin Body Lotion, MoliCare® Hand Cream, MoliCare® Protection Foam

# 3 easy steps to select the right pad



## STYLE

Various styles suit different levels of mobility.

### Active


Most styles suitable. Selection based on personal preference.



MoliCare Premium Pads and Pants

### Active with restricted mobility and/or cognitive difficulties (e.g. Dementia)

Pull-up recommended. Looks and feels like normal underwear.



MoliCare Premium Mobile

### Restricted/supported mobility or bedbound

All-in-one product recommended. Easiest option for the carer.



MoliCare Premium Elastic, Premium Form and Fixpants

## ABSORBENCY

Giving the right level of protection — Choose an absorbency that suits your level of incontinence and keeps you dry.



Less Absorbent



More Absorbent

Look for the drop rating on pack  
**More drops = Higher absorbency**

## SIZE

To ensure maximum comfort and leakage protection:

- Measure waist at the largest width between the waist and hips.
- Select the product size based on the size guide printed on the product packaging.



Leakage may occur if the product is too large or too small.





# MoliCare® Assortment Overview


	MoliCare Premium Lady Pad	MoliCare Premium Men Pad	MoliCare Pad	MoliCare Rectangular	MoliCare Premium Mobile	MoliCare Premium Form	MoliCare Premium Elastic	MoliCare Slip XS	MoliCare Premium Fixpants Long Leg	MoliCare Premium Fixpants Short Leg
										
Gender and Mobility Level										
Absorbency Level										
1										
2										
3										
4										
5										
6										
7										
8										
9										
10										
Special range										


Due to different cuts, shapes and sizes of the absorbent core in each product group, the exact absorption level differs within each drop level, (i.e. 1 drop does not equal a defined ISO absorbency). For the exact absorbency measures, refer to specific product information throughout the catalogue. Colours correspond to colour coding of packaging.


Category symbols


  
Unisex products


  
Products specifically for female wearers


  
Products specifically for male wearers


  
Available sizes

  
Available absorbency levels

  
Products specifically for active wearers

  
Products for wearers with restricted mobility


  
Products for wearers with supported mobility

  
Products for bed bound wearers



# MoliCare® Fitting Guide

## Find the right fit



Due to different cuts, shapes and fixation systems, different product categories follow different sizing guides. Measure waist at the largest width between the waist and hips and select the product size based on the size guide printed on the product packaging or below.

Product Group	Absorbency	cm Size	35	40	45	50	55	60	65	70	75	80	85	90	95	100	105	110	115	120	125	130	135	140	145	150	155	160	165	170	175	180	185	190	195	200	200+
MoliCare Premium Mobile	6	XS			XS (45 - 70cm)																																
	5 6 8	S						S (60 - 90cm)																													
	5 6 8 10	M										M (80 - 120cm)																									
		L														L (100 - 150cm)																					
	6 8 10	XL																			XL (130 - 170cm)																
MoliCare Premium Elastic	6 8 9 10	S								S (70 - 90cm)																											
		M										M (85 - 120cm)																									
		L															L (115 - 145cm)																				
		XL																			XL (140 - 175cm)																
MoliCare Slip XS	6	XS		XS (40 - 70cm)																																	

MoliCare Premium Fixpants		cm Size	35	40	45	50	55	60	65	70	75	80	85	90	95	100	105	110	115	120	125	130	135	140	145	150	155	160	165	170	175	180	185	190	195	200	200+
Long Leg	Short Leg																																				
-	Yellow	S			S (45 - 80cm)																																
Blue		M						M (60 - 100cm)																													
Brown		L										L (80 - 120cm)																									
Green		XL														XL (100 - 160cm)																					
Grey	-	XXL																					XXL (140 - 180cm)														
Orange		XXXL																										XXXL (160 - 200cm)									
Purple		5XL																																			5XL (200 - 240cm)



# MoliCare Advanced Features



## Active Skin Protection

Skin-friendly, pH balance ensures maximum skin compatibility. Keeps skin dry with a 3-layer absorbent core and advanced fast dry effect.



## Odour neutraliser

Neutralises unpleasant odours for enhanced well-being.



## Dermatologically tested

Helps to maintain healthy skin.



## Adjusted fit

Safety & comfort through optimal fit; elasticised side panels ensure better adaptation to body shape and movements, delivering improved comfort.



## Air active

Skin-friendly due to high quality air active material with an extremely low potential to trigger allergic reactions.



## Anatomically shaped fit

MoliCare Premium Men Pants: Feels like regular underwear. Front shifted absorbent core for extra protection where men need it most.



## Anti-leakage system

For an improved feeling of safety.



## Comfortable fit

The soft waistband, the skin-friendly, air active material and the cross-elasticised threads ensure high wearing comfort.



## Discreet underwear feeling

Discreet underwear-feel through anatomical shape and soft, stretchy material for a comfortable fit.



## Elasticated side panels

Elasticised side panels in the hip area ensure better adaptation to body shape and movements, delivering improved comfort for the patients.



## Ergonomic use

Approximately 20% faster application as well as easier and more ergonomic handling for caregivers.<sup>1</sup>



## Extra protection

MoliCare Premium Bed Mat: Absorbs liquid and provides extra protection for various surfaces.



## Hygienic care

Single-use secure protection.



## Prefolded & compressed

For easy use and storage.



## Re-closable tapes

2 re-closable tapes (one for each side) for easy and secure closing.



## Secure fixation

Designed to be worn with regular underwear.



## Secure fixation

MoliCare Premium Fixpants: Ensure a safe and hygienic fixation of anatomically shaped large incontinence pads.



## Soft & skin-friendly

Soft, breathable material for healthy skin.



## Underwear like

Ideal body fit with cross-elasticated threads and extra long leg joints to ensure high wearing comfort.



## Washable

Washable up to the number of times indicated on the label in a laundry net and can be tumble dried at low heat.



## Wetness indicator on the back sheet

Shows when it is time to change.



## 3-layer absorbent core

Improved absorption core technology that helps to prevent skin irritations.



## Advanced skin protection

Thanks to a balanced pH.



## Anti-leakage system

For an improved feeling of safety.



## Dermatologically tested

Helps to maintain healthy skin.



## Extra odour neutralisation

Superabsorbent polymers neutralise unpleasant odours



## Active skin protection

Up to 64% faster urine absorption compared to the previous MoliCare® Premium Form.

<sup>1</sup> W. Schneider, H. Mutschler, R. Haller, C. le Mathe: Test for ergonomic evaluation of incontinence products, Oct 5-8, 2015 Dijon. The reduction from 4 to 2 tapes (one for each side) combined with no side panels in the front simplifies the application by reducing the steps from 8 to 4 compared to a traditional slip with 4 side panels and 4 tapes. Additionally, due to the elasticated side panels, there is less need of readjusting the tapes.







# PADS

1 – 5 

## MoliCare® Premium Lady Pad



1 Drop	USL Code	HARTMANN Code	Fit	ISO Max Absorbency in mL*	Units Per Pack	Pack Per Carton
	16479	168132	Women	181	14	12
2 Drops	USL Code	HARTMANN Code	Fit	ISO Max Absorbency in mL*	Units Per Pack	Pack Per Carton
	8959	168634	Women	331	14	18
3 Drops	USL Code	HARTMANN Code	Fit	ISO Max Absorbency in mL*	Units Per Pack	Pack Per Carton
	8960	168644	Women	504	14	12
5 Drops	USL Code	HARTMANN Code	Fit	ISO Max Absorbency in mL*	Units Per Pack	Pack Per Carton
	1232	168670	Women	1029	14	12

\* Measurement in accordance with ISO 11948/1 (Rothwell method)



MoliCare® Premium Lady Pad 1 Drop




MoliCare® Premium Lady Pad 2, 3, 5 Drops




# PADS




## MoliCare® Premium Men Pad



2 Drops	USL Code	HARTMANN Code	Fit	ISO Max Absorbency in mL	Units Per Pack	Pack Per Carton
	1381	168600	Men	330	14	12

4 Drops	USL Code	HARTMANN Code	Fit	ISO Max Absorbency in mL	Units Per Pack	Pack Per Carton
	1383	168705	Men	546	14	12



## PADS

2 – 4 

### MoliCare® Pad



2

Drops



USL  
Code

1208

HARTMANN  
Code

168100

Fit

Unisex

ISO Max  
Absorbency in mL

321

Units  
Per Pack

28

Pack  
Per Carton

12

3

Drops



1209

168101

Unisex

481

28

6

4

Drops



1210

168102

Unisex

865

28

6

## SUPPORT PADS

3 

### MoliCare® Rectangular



3

Drops



USL  
Code

7165

HARTMANN  
Code

160890

Dimensions

20 x 40 cm

ISO Max  
Absorbency in mL

430

Units  
Per Pack

50

Pack  
Per Carton

6



# PULL-UPS

5 - 10 XS - XL

## MoliCare® Premium Mobile



DISCREET & COMFORTABLE



### 5 Drops



USL Code	HARTMANN Code	Size	Waist Measurement	ISO Max Absorbency in mL	Units Per Pack	Pack Per Carton
9247	915861	S	60 - 90 cm	941	14	4
9248	915852	M	80 - 120 cm	1017	14	3
9249	915853	L	100 - 150 cm	1198	14	4

### 6 Drops



9250	915840	XS	45 - 70 cm	1361	14	4
1226A	915831	S	60 - 90 cm	1475	14	4
1227A	915832	M	80 - 120 cm	1662	14	3
1228A	915833	L	100 - 150 cm	1963	14	4
9251	915834	XL	130 - 170 cm	2140	14	4

### 8 Drops



9252	915871	S	60 - 90 cm	1791	14	4
9253	915872	M	80 - 120 cm	2015	14	3
9254	915873	L	100 - 150 cm	2279	14	4
92550	915874	XL	130 - 170 cm	2421	14	4

### 10 Drops



16260	915878	M	80 - 120 cm	2626	14	3
16261	915879	L	100 - 150 cm	2613	14	4
16262	915880	XL	130 - 170 cm	2757	14	4

# LARGE PADS


4 – 9 

## MoliCare® Premium Form



Recommended with MoliCare® Premium Fixpants 

4 Drops	USL Code	HARTMANN Code	Fit	ISO Max Absorbency in mL	Units Per Pack	Pack Per Carton
	15622	168404	Unisex	1493	32	4

5 Drops	USL Code	HARTMANN Code	Fit	ISO Max Absorbency in mL	Units Per Pack	Pack Per Carton
	15624	168405	Unisex	1662	32	4

6 Drops	USL Code	HARTMANN Code	Fit	ISO Max Absorbency in mL	Units Per Pack	Pack Per Carton
	15626	168406	Unisex	2353	32	4

6 Drops 	USL Code	HARTMANN Code	Fit	ISO Max Absorbency in mL	Units Per Pack	Pack Per Carton
	15634	168401	Men	2441	32	4

7 Drops	USL Code	HARTMANN Code	Fit	ISO Max Absorbency in mL	Units Per Pack	Pack Per Carton
	15628	168407	Unisex	2598	32	4

8 Drops	USL Code	HARTMANN Code	Fit	ISO Max Absorbency in mL	Units Per Pack	Pack Per Carton
	15630	168408	Unisex	3017	32	4

9 Drops	USL Code	HARTMANN Code	Fit	ISO Max Absorbency in mL	Units Per Pack	Pack Per Carton
	15632	168409	Unisex	3497	16	4



## LARGE PADS

8 +SIZE 

### MoliCare® Premium Form + SIZE



Recommended  
with  
MoliCare®  
Premium Fixpants 

8 Drops	USL Code	HARTMANN Code	Fit	ISO Max Absorbency in mL	Units Per Pack	Pack Per Carton	8
	1212	168519	Unisex	3347	18	3	

### MoliCare® Premium Form STOOL



Recommended  
with  
MoliCare®  
Premium Fixpants 

STOOL	USL Code	HARTMANN Code	Fit	ISO Max Absorbency in mL	Units Per Pack	Pack Per Carton
	1213	168419	Unisex	1295	32	4

## FIXPANTS LONG LEG

M - 5XL 

### MoliCare® Premium Fixpants Long Leg



Recommended with MoliCare® Premium Form 



USL Code	HARTMANN Code	Size	Waist Measurement	Colour	Units Per Pack	Pack Per Carton
14528	947796	M	60 - 100 cm	Blue	5	40
1372A	947791	M	60 - 100 cm	Blue	25	8
14529	947797	L	80 - 120 cm	Brown	5	40
1378A	947792	L	80 - 120 cm	Brown	25	8
14530	947798	XL	100 - 160 cm	Green	5	40
1382A	947793	XL	100 - 160 cm	Green	25	8
14536	947799	XXL	140 - 180 cm	Grey	5	40
1379A	947794	XXL	140 - 180 cm	Grey	25	8
14537	947815	XXXL	160 - 200 cm	Orange	5	40
1371A	947789	XXXL	160 - 200 cm	Orange	25	8
15636	947830	5XL	200 - 240 cm	Purple	5	25
1377	947825	5XL	200 - 240 cm	Purple	25	5



## FIXPANTS SHORT LEG

S - XL 

### MoliCare® Premium Fixpants Short Leg



Recommended with MoliCare® Premium Form 



USL Code	HARTMANN Code	Size	Waist Measurement	Colour	Units Per Pack	Pack Per Carton
14538	947759	S	45 - 80 cm	Yellow	5	60
1373A	947711	S	45 - 80 cm	Yellow	25	12
14539	947760	M	60 - 100 cm	Blue	5	60
1374A	947712	M	60 - 100 cm	Blue	25	12
14560	947761	L	80 - 120 cm	Brown	5	60
1375A	947713	L	80 - 120 cm	Brown	25	12
14566	947762	XL	100 - 160 cm	Green	5	60
1376A	947714	XL	100 - 160 cm	Green	25	12

# ALL-IN-ONE SLIPS


6 - 10 S - XL 

## MoliCare® Premium Elastic




COMFORTABLE  
& EASY USE



6 Drops	USL Code	HARTMANN Code	Size	Waist Measurement	ISO Max Absorbency in mL	Units Per Pack	Pack Per Carton
	9125	165271	S	70 - 90 cm	1667	30	3
	9126	165272	M	85 - 120 cm	2159	30	3
	9127	165273	L	115 - 145 cm	2349	30	3
	7548	165274	XL	140 - 175 cm	2786	14	4

8 Drops	USL Code	HARTMANN Code	Size	Waist Measurement	ISO Max Absorbency in mL	Units Per Pack	Pack Per Carton
	9128A	165471	S	70 - 90 cm	2356	26	3
	9129A	165472	M	85 - 120 cm	3144	26	3
	9130A	165473	L	115 - 145 cm	3220	24	3
	14068	165474	XL	140 - 175 cm	3591	14	4

9 Drops	USL Code	HARTMANN Code	Size	Waist Measurement	ISO Max Absorbency in mL	Units Per Pack	Pack Per Carton
	14411A	165571	S	70 - 90 cm	2546	26	3
	14412A	165572	M	85 - 120 cm	3521	26	3
	14413A	165573	L	115 - 145 cm	3856	24	3
	1380	165574	XL	140 - 175 cm	3951	14	4


10 Drops	USL Code	HARTMANN Code	Size	Waist Measurement	ISO Max Absorbency in mL	Units Per Pack	Pack Per Carton
	14680A	165671	S	70 - 90 cm	2771	22	3
	14681	165672	M	85 - 120 cm	3696	14	4
	14682	165673	L	115 - 145 cm	4279	14	4
	14683	165674	XL	140 - 175 cm	4380	14	4




## ALL-IN-ONE SLIPS

6 XS 

### MoliCare® Slip XS




6 Drops	USL Code	HARTMANN Code	Size	Waist measurement	ISO max absorbency in mL	Units per pack	Pack per carton
	15640	169190	XS	40 - 70 cm	1367	30	4

## BED MATS

5 – 9 

### MoliCare® Bed Mat Eco





5 Drops	USL Code	HARTMANN Code	Dimensions	ISO Max Absorbency in mL	Units Per Pack	Pack Per Carton
	1205	159001	40 x 60 cm	497	300	1


### MoliCare® Premium Bed Mat



7 Drops	USL Code	HARTMANN Code	Dimensions	ISO Max Absorbency in mL	Units Per Pack	Pack Per Carton
	1206	161066	40 x 60 cm	561	150	1
	1207	161070	60 x 90 cm	1241	100	1

8 Drops	USL Code	HARTMANN Code	Dimensions	ISO Max Absorbency in mL	Units Per Pack	Pack Per Carton
	8204	1610880	60 x 90 cm	1407	30	4

9 Drops	USL Code	HARTMANN Code	Dimensions	ISO Max Absorbency in mL	Units Per Pack	Pack Per Carton
	8203A	161074	40 x 60 cm	1168	30	3

 40 x 60 cm variant not compressed



Helps. Cares. Protects.

PAUL HARTMANN Pty Ltd. ABN 35 000 099 589.  
Level 5, 1 Thomas Holt Drive, Macquarie Park NSW 2113  
PO Box 943, North Ryde BC NSW 1670  
® Registered trademark PAUL HARTMANN AG.

9166444\_052023



For further information on MoliCare continence management products, please contact your USL Account Manager or call customer service on 0800 658 814.