

PRODUCT SPECIFICATION

Product:	Nitrile Disposable Glove
Article No.:	WRP NBR 3549
Style no:	WRP NBR 3549 - E1 35
Material:	Nitrile Rubber
Colour:	Blue, Violet Blue, White
Size:	XS, S, M, L, XL, XXL
Texture:	Finger tip textured
Steriisation:	Non sterile
Applications:	Food handling, catering, cleaning and janitorial, light industrial chemical protection, agriculture, printing, construction, automotive industries, electroplating, environmental waste clean up, house hold chores, medical and healthcare

PHYSICAL DIMENSION:

Thickness at palm : mm	0.05 (min)
-------------------------------	------------

Width at palm (mm):	XSmall: 80 +/- 10
	Small: 85 +/- 10
	Medium: 95 +/- 10
	Large: 110 +/- 10
	Xlarge: 115 +/- 10
	XXLarge: 125 +/- 10

Length (mm):	240 +/- 10
---------------------	------------

Weight / piece (g):

Xsmall:	2.5 +/- 0.3
Small:	3.0 +/- 0.3
Medium:	3.5 +/- 0.3
Large:	4.0 +/- 0.3
Xlarge:	4.5 +/- 0.3
XXLarge:	5.0 +/- 0.3

TECHNICAL SPECIFICATION

Physical Properties	Unaged	Aged (22hr, 100c)
Tensile Strength at break	14MPa (min)	14 MPa (min)
Elongation at break	500%	500%

COMPLIANCE & TESTED IN ACCORDANCE TO	ASTM : D 6319- 10 (2015) , EN 455 : Part 1,2,3,4 , EU CE Regulation 2016/ 425, FDA 510 K, EN 420 : 2003 + A1 2009 , EN ISO 374-1: 2016, EN 374-2: 2014, EN 374-4: 2013, EN ISO 374-5: 2016, EN 16523: 2015, EU 10/2011
---	--

MANUFACTURED BY: WRP ASIA PACIFIC Sdn Bhd

ADDRESS: Lot 1, Jalan 3, Kawasan Perusahaan Bandar Baru Salak Tinggi,
43900 Sepang,
Selangor Darul Ehsan,
Malaysia

TELEPHONE: + 6 0687061486

FAX + 6 0687061557

EMAIL ADDRESS customer_wrp@wrpworld.com; sales@wrpworld.com

WEBSITE www.wrpworld.com