



Infected wounds?

Treat with **Atrauman® Ag**



**Impregnated Silver Wound Contact Layer from HARTMANN.**

Atrauman® Ag is a low adherent wound contact layer which contains silver. Designed specifically for contaminated, traumatic wounds helping reduce the risk of infection.

**Reduces the Risk of Bacterial Burden.**

Atrauman® Ag is coated in silver which has been proven to minimise the bacterial burden on wounds. The silver ions bind with bacteria and destroy them within the dressing whilst protecting the wound bed.

# For a wide range of indications

## PRODUCT BENEFITS

**Effectively kills** a broad spectrum of bacteria

**Lowers toxicity** to human cells thanks to reduced release of silver ions <sup>[1]</sup>

**Fine weave** allows passage of exudate away from the wound surface

**Improved wound** healing outcomes

**Infection prevention** in wound treatment

## User considerations

### Suggested Wear Time

Up to 3 days


### Suggested secondary dressings

- **Biatain® / Dry Max Foam** - for moderate exuding wounds
- **Zetuvit® Plus / Dry Max Extra Soft** - for highly exuding wounds

### Suggested Fixation

- **Peha-haft®** - cohesive bandage
- **Omnifix®** - non woven tape
- **Hydrofilm®** - waterproof film

### Application



Skin Tears



Infected Post Op Wounds




Grazes



Diabetic Foot Ulcers



Arterial Ulcers



Venous Leg Ulcers

CLICK  
HERE to  
link back  
to website  
listing



PRODUCT	SIZE	UNITS PER BOX	CODE
Atrauman®Ag sterile, individually sealed	5 x 5 cm	10pcs	8904
	10 x 10 cm	10pcs	8906
	10 x 20 cm	10pcs	8908

<sup>[1]</sup> K. Ziegler et al. (2006). Reduced Cellular Toxicity of a New Silver Containing Antimicrobial Dressing and Clinical Performance in Non-Healing Wounds. *Skin Pharmacol Physiol* 2006;19:140-146