

The **Dyna-Flex Cushion** incorporates the most exclusive combination of visco elastic & CME foam technology, providing not only outstanding pressure re-distribution, but great patient comfort. The temperature sensitive memory foam softens at patient interface, thus giving even weight distribution, whilst the CME underside prevents the foam from bottoming out. Suitable for patients considered to be at **"High Risk"** of pressure ulcer development, **Dyna-Flex** gives nursing security without the patient challenging the issue of comfort. Complete with a well designed cover and non-slip base, this product is ideal for use in all hospital, nursing or homecare environments.

	Features	Benefits	Technical Specific	fications	
1	Temperature sensitive memory foam	excellent comfort & pressure re-distribution properties Risk Category	Product Code Risk Category	CUSH/FLEX/43/43/7.5 High Risk	
2	Multi-stretch, vapour permeable fabric (impregnated with an anti fungal agent)	Reduces shear & friction forces, whilst assisting infection control policies	Weight Limit Warranty Dimensions Foam Density Fatigue Class	24 Stone / 152kg 2 Years (Length) 43cm x (Width) 43cm x	
3	High frequency welded cover	Prevents fluid ingress and aids decontamination processes (can be laundered up to 95 °C)		(Height) 7.5cm (17" x 17" x 3") Variable Sizes Available Upon Request Visco – 55–60kg/m³	
4	Non-slip base	Prevents slippage, encouraging patient safety		Nominal Hardness – Soft to Medium CME Base – 38–40kg/m³ Nominal Hardness – 200 Newtons	
(5)	Fully replaceable parts	Ensures great value for money		Severe	
6	Available with sag compensator	Suitable for wheelchair users	Fire Retardancy	BS7176; 2007	
			Weight	2kg	
			Colour	Wine	
			Hypochlorite Solution compatible	1,000 ppm diluted or 0.1% (or to be used in conjunction with local infection control guidelines)	