



Hillrom™


Welch Allyn® MacroView® Otosopes

EXPAND YOUR VIEW

Experience bigger, brighter views with our Welch Allyn MacroView and MacroView Plus Otoscopes. Enhance exams with the MacroView Plus Otoscope, and the option to capture and share ear images for consultation. Both devices work on existing Welch Allyn 3.5 V power sources.

Pair the MacroView otoscope with Welch Allyn LumiView™ Clear Single-Use Ear Specula for brighter views.

Save and share ear images with the MacroView Plus device when combined with our iExaminer® SmartBracket™ and iExaminer Pro App.




| Specifications | MacroView Plus Otoscope for iExaminer | MacroView Basic Otoscope |
|---------------------------------------|--|--------------------------|
| Model Number | 238-3 | 238-2 |
| iExaminer digital-ready | Yes | No |
| Warranty | 5 Years | 1 Year |
| Height x Width x Depth | 7 cm h x 4.5 cm w x 5.5 cm d | |
| Weight | 115 grams | |
| Sealed Pneumatic otoscopy | Yes | |
| Illumination type | Ring fiber optics | |
| Lamp Type | SureColor™ LED | |
| Lamp Life, rated | 60+ years based on typical use ¹ | |
| Field of view | ~3 times that of a single-lens otoscope, like the standard-style 250-2. | |
| Magnification | ~5X | |
| Focus-free design | Extended focal length virtually eliminates the need to focus on wide-view otoscope | |
| Total light output | 5200 lux; 2.5 times brighter than traditional halogen scopes. | |
| Color | Color rendering index (CRI): 95 typical; R9: 90 typical; Color Temperature: 3400K - 4500K | |
| LED warranty | 7 years | |
| Optimized for the 719-3 Li-Ion handle | Stays bright over the life of a battery charge (unlike older designs that dim with use); no color change with dimming (unlike halogen) | |

3X LARGER Viewing Area



Standard
Otoscope



MacroView
Otoscope

| Specifications | Welch Allyn® Lithium Ion Plus Power Handle with USB Charger | Welch Allyn® Universal Desk Charger |
|-------------------------------|---|---|
| Model Number | 719-3 | 719-DSK |
| Height x Width X Depth | 14 cm x 3 cm x 3 cm | 10 cm x 19 cm x 11 cm |
| Weight | 140 grams | 1.2 kg |
| Compatibility | 719-DSK Welch Allyn® Universal Desk Charger; 71140 chargers | 71140, 719-3, 71910, 71670 handles. 12800 & 22800 used with 71249 Adapter Sleeve |
| Integrated Charging | USB-C | — |
| Warranty | Two years | One year |
| On-time | 5.5 hours; Approximately 660 30-second exams | — |
| Battery Type | Lithium Ion | — |
| USB Quick-charge time (0-75%) | 90 minutes | — |
| Power Cord and Cord Options | 719-USB | #18 AWG, 3-wire grounded cord, 6 feet long medical grade detachable with IEC60320 inlet |
| Input | — | 100-240V, 50-60 Hz, 100mA max, Class I, Continuous Operation |
| Output | — | 2.4V - 4.0V DC 80mA max |

SELECT MACROVIEW ACCESSORIES

52432-CLR-1: 2.75 mm LumiView Single-Use Clear Ear Specula, Pediatric (case)

52434-CLR-1: 4.25 mm LumiView Single-Use Clear Ear Specula, Adult (case)

11840-TST-US: iExaminer Resolution Test Card

11840-IEX-US: iExaminer SmartBracket Accessory

21504: Insufflation Bulb

24320: Welch Allyn® SofSeal™ Specula Tip (supports insufflation)



Distributed in New Zealand by:

0800 658 814 | uslcustomerservices@uslmedical.co.nz | uslequipment.co.nz

