Suction Tubing

Australian Made, high quality and cost effective

Fairmont Medical has excelled in providing reliability of supply to meet our customers needs, providing an extensive range of suction tubing solutions.

Red Connector


- 100% Australian Made
- Trumpeted opening allows for variety of suction devices to be attached
- Smooth inner bore surface allows for easy flow of suction

Sterile

Product Code	Length	Tube ID	Tube OD	Box Quantity
DST3002	2m	6mm	9mm	100
DST3003	3m	6mm	9mm	60
DST3045	4.5m	6mm	9mm	40

Non-Sterile

Product Code	Length	Tube ID	Tube OD	Box Quantity
DST3012	2m	6mm	9mm	100
DST3013	3m	6mm	9mm	60
DST3145	4.5m	6mm	9mm	40



White Hexagonal Connector

- Universal Connector
- Provides quick and easy connection without wastage

Sterile

Product Code	Length	Tube ID	Tube OD	Box Quantity
DST3202	2m	6mm	9mm	100
DST3203	3m	6mm	9mm	60
DST3245	4.5m	6mm	9mm	40
DST3206	6m	6mm	9mm	40
DST3208	8m	6mm	9mm	25

Non-Sterile

Product Code	Length	Tube ID	Tube OD	Box Quantity
DST3212	2m	6mm	9mm	100
DST3213	3m	6mm	9mm	60
DST3214	4.5m	6mm	9mm	40
DST3215	1.5m	6mm	9mm	100
DST3216	6m	6mm	9mm	40
DST3218	8m	6mm	9mm	25



Specialty Suction Tubing

Nelson ENT Suction Tubing

• Tubing with Red Connector and silicone tubing

Sterile

Product Code	Length	Tube ID	Tube OD	Silicone Tubing	Box Quantity
DST3030	3m	6mm	9mm	200mm	60
DST3035	5m	6mm	9mm	500mm	40



Soft Tubing

1 White Hexagonal Connector

Sterile

Product Code	Length	Tube ID	Tube OD	Box Quantity
DST3302	2m	6mm	9mm	45
DST3303	3m	6mm	9mm	50
DST3304	4m	6mm	9mm	45




Soft Tubing

• 1 White Hexagonal Connector with Double Barbed Connector

Sterile

Product Code	Length	Tube ID	Tube OD	Box Quantity
DST3310	1m	5.5mm	8.8mm	100




Soft Tubing

• 2 White Hexagonal Connectors

Sterile

Product Code	Length	Tube ID	Tube OD	Box Quantity
DST3305	4m	6mm	9mm	45
DST3345	4.5m	6mm	9mm	45




Specialty Suction Tubing

Soft Tubing

No Connectors

Sterile

Product Code	Length	Tube ID	Tube OD	Box Quantity
DST3300	1.5m	6mm	9mm	100



Wide Bore Tubing

• 2 White Hexagonal Connectors

Sterile


Product Code	Length	Tube ID	Tube OD	Box Quantity
DST3402	2m	6.9mm	10.9mm	60
DST3403	3m	6.9mm	10.9mm	60
DST3445	4.5m	6.9mm	10.9mm	40
DST3406	6m	6.9mm	10.9mm	30
DST3408	8m	6.9mm	10.9mm	25



Irrigation Tubing

Sterile

Product Code	Length	Tube ID	Tube OD	Box Quantity
DIT5015	150mm	5mm	7.6mm	50



Blue Frazier Connector

• 1 Blue Frazier Connector and 1 White Hexagonal Connector

Sterile

Product Code	Length	Tube ID	Tube OD	Box Quantity
DST3032	2m	6mm	8.6mm	80
DST3303	3m	6mm	8.6mm	60
DST3034	4m	6mm	8.6mm	40
DST3036	6m	6mm	8.6mm	40
DST3038	8m	6mm	8.6mm	25



Specialty Suction Tubing

Caws Tubing

Double/Universal Connector

Sterile

Product Code	Length	Tube ID	Tube OD	Box Quantity
DST3610	100mm	6mm	9mm	100
DST3615	500mm	6mm	9mm	50



Double/Universal Connector

Sterile

Product Code	Length	Tube ID	Tube OD	Box Quantity
DST3720	200mm	4.9mm	6.9mm	20

• Silicone tubing with double/universal connector

Sterile

Product Code	Length	Tube ID	Tube OD	Box Quantity
DST3730	300mm	4.8mm	9.5mm	50
DST3750	500mm	4.8mm	9.5mm	50



Lidia Tubing

• With Male Luer Connector

Sterile

Product Code	Length	Tube ID	Tube OD	Box Quantity
DST3603	3m	5.9mm	8.9mm	60



Wide Bore Liposuction Tubing

Blue Connector

Sterile

Product Code	Length	Tube ID	Tube OD	Box Quantity
DST3503	3m	9.8mm	15mm	10

