

- Sterile, double wrapped in individual peel pouches
- Available with or without VC (Vacuum control) vent
- Available with or without sump tip to reduce the risk of tissue damage

## **Derby Yankauers**

Clear, acrylic suction devices, similar to traditional yankauer suckers. With or without sump tip to reduce potential tissue suction during use.

Product Code	Туре	Length	CH	Inside Diam. mm	Qty. Per Box
YS-3004	with sump tip	265mm	24CH	8	50
YS-3005	with sump tip and VC	265mm	24CH	8	50
YS-3014	straight tip	265mm	24CH	8	50
YS-3015	straight tip and VC	265mm	24CH	8	50
YS-3000	Midi Yankauer with straight tip	265mm	20CH	6.66	50
YS-3001	Midi Yankauer with VC	265mm	20CH	6.66	50

## **Derby Yankauers Sets**

Derby Yankauer Suckers pre-attached to 2m or 3m connecting tube. Available with or without 'sump' tip to reduce the risk of tissue damage during use.

- All Pennine suction sets come sterile, double wrapped in individual peel pouches
- Available with or without VC (vacuum control) hole

Set with 2m tubing	CH/Tip	Qty. Per Box	Set with 3m tubing	Qty. Per Box	Tube Diam. mm	Yankauer type
SYS-5820	24CH	20	SYS-5830	10	6	YS-3004
SYS-5821	24CH	20	SYS-5831	10	6	YS-3005
SYS-5520	24CH	20	SYS-5530	15	6	YS-3014
SYS-5521	24CH	20	SYS-5531	15	6	YS-3015



## -

- Switch system allows the operator to manually open and close suction
- Easy connection to Pennine range of tubing

## Selectavac Yankauer

Product Code	Туре	Length	CH/Tip	Inside Diam. mm	Qty. Per Box
YS-3009	Bulb tip	265	24CH	8	50
YS-3010	Straight tip	265	24CH	8	50