



Lastolan: the strong, reusable short-stretch bandage for oedema reduction and support.





A very strong cotton crepe compression bandage.



Lastolan provides high compression for conditions in which venous flow is compromised.



What is Lastolan?

Lastolan is a very strong, short-stretch, 100% cotton crepe bandage with approximately 90% extensibility. Lastolan is extraordinarily durable and generates a high working pressure support and a low resting pressure. This allows decongestion of even deep sub-fascial venous areas.

Lastolan adapts to changes in circumference caused by oedema subsidence, so can remain in place for 1–2 days, except in the acute phase of severe oedema.

Lastolan is air-permeable and kind to the skin. It can be washed at temperatures up to 85°C, and steam sterilised. Its durability allows for multiple re-use.

What is Lastolan used for?

Lastolan is ideal for treating conditions that require a high working pressure support, including:

- lymphoedema
- acute and chronic venous stasis
- venous leg ulcers
- thrombophlebitis
- post-sclerotherapy
- musculoskeletal injuries
- for most common compression techniques such as Pütter, Fischer and Sigg

For moderate compression, application of a 2-layer Lastolan bandage allows compression of deep sub-fascial veins.

	Width	Unit	Code No	Units per case
Lastolan stretched length 5m. Individually boxed.	6cm	10	931 881	32
	8cm	10	931 882	24
	10cm	10	931 883	18
	12cm	10	931 884	16
Lastolan stretched length 5m. Box of 10 bandages.	6cm	10	931 891	50
	8cm	10	931 892	36
	10cm	10	931 893	30
	12cm	10	931 894	24

Lastolan: the benefits.

- Attains high working pressure
- Very low resting pressure
- Air permeable
- Durable – multiple use possible
- Reusable and economical due to washability
- Low slip when applied with Pütter technique

THE HARTMANN ALLIANCE

NEW SOUTH WALES

Paul Hartmann Pty Ltd
Phone: 1800 805 839

VICTORIA & TASMANIA

Sutherland Medical Pty Ltd
Phone: 1300 664 027

QUEENSLAND

Medical & Surgical Requisites Pty Ltd
Phone: 1800 801 226

SOUTH AUSTRALIA & NT

McNeil's Surgical Pty Ltd
Phone: 08 8363 0888

WESTERN AUSTRALIA

Perth Surgical Supply Co. Pty Ltd
Phone: 08 9344 3111

HARTMANN

